
CONGEN

SureTaq Hotstart Taq Polymerase II

Art. No. F4003
100 rxn

User Manual

September 2025

 Inhalt

1 Allgemeines3

1.1 Beschreibung.....3

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1 General Information

1.1 Description

SureTaq Hotstart Taq Polymerase II is an additional component for use with SureFast® / SureFood® real-time kits.

SureTaq Hotstart Taq Polymerase II is not verified lot-specifically with SureFast® / SureFood® real-time Kits.

1.2 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
T	Taq Polymerase	1 x 80 µl	Dark Red

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.3 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP products / SureFood® PREP products)
- Real-time PCR methods (e.g. SureFast® products / SureFood® products)

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Further information about the calculation of the total number of reactions needed (samples and control reactions) for the specific PCR assay is given in the product-related documents of the SureFast® / SureFood® real-time PCR products.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

3 Further Information

3.1 Technical Support

For further questions please send an e-mail to info@congen.de.