CONGEN

SureFast[®] Master Probe

Art. No. F4006-2
500 rxn

User Manual

March 2021

 **Inhalt**

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1 General Information

1.1 Description

The SureFast® Master Probe can be used for the specific detection of PCR products based on probe technology. The kit includes all required reagents except of primer, probe and nucleic acid.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, e.g. Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, R-Biopharm RIDA®CYCLER and Agilent AriaMx.

1.2 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	5 x 1300 µl	Yellow
2	PCR Water	5 x 450 µl	White
3	Taq Polymerase	5 x 11 µl	Red

1.3 Additionally required equipment and materials

- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- specific primer, probe, DNA
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.4 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control. Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
Primer 1* (not included in the kit)	1.7 µl	18.7 µl
Primer 2* (not included in the kit)	1.7 µl	18.7 µl
Probe* (not included in the kit)	0.6 µl	6.6 µl
PCR Water	3.4 µl	37.4 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

*Recommendation: concentration: 5 pmol/µl

**Duplex- or Multiplex-PCR systems maybe need a higher concentration of Taq Polymerase. For Duplex- or Multiplex-PCR reactions we recommend the use of double amount (0.2 µl/rxn). In this case the SureFast® Master Probe kit suffices for 250 reactions. Extra Taq Polymerase can be order directly with Art. No. F4004 / F4005 SureTaq.

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions have to show the correct results.

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.