**SureFast® Master Probe PLUS (100 Reakt.)**

Art. Nr. F4007-2 (5\*100 Rkt.)

Januar 2017

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**Description**

The SureMaster Probe PLUS Kit can be used for the specific detection of PCR products based on probe technology. Each reaction contains an internal amplification control (PLUS). The kit is optimized for detection of PCR products in a final volume of 25 µl. The reaction mix includes all necessary reagents except to the specific primer, probe (recommendation reporter FAM, quencher TAMRA or dark-quencher) and DNA for the detection system. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions at 522 nm and 553 nm at the same time (e.g. Roche LightCycler® 480, Rotor-Gene Q, BioRad CFX96, Agilent MxSeries etc.).

**Limit of Detection**

The internal amplification control system is a synthetic one. There are no comparable natural sequences known. It is possible that parts of the sequences used in the detection system are conforming to the sequence of the internal control system. In this case the detection limit of the detection system will be influenced. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

**Kit components and storage**

Kit Code	Reagent	Amount	Lid Color
1	Probe Mix	1 x 1400 µl	Yellow
2	PCR Water	1 x 950 µl	White
3	Taq Polymerase	1 x 11 µl	Red

Store all reagents at -20°C and protected from light.

**Additionally required equipment and materials**

- Real-time PCR instrument with two detection channels (522 nm, 553 nm)
- Real-time PCR consumables (plates, tubes, foils)
- Pipettes with filter tips
- Specific primer, probe, DNA

**Protocol**

## 1. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix contains an internal amplification control per reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
SureMaster Probe PLUS Reaction Mix	12.5 µl	137.5 µl
Primer 1*	1.7 µl	18.7 µl
Primer 2*	1.7 µl	18.7 µl
Probe*	0.6 µl	6.6 µl
PCR Water	3.4 µl	37.4 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>
Volume DNA	5 µl	

\*Recommendation concentration: 5 pmol/µl

Mix each master-mix well and centrifuge shortly before use.

2. Setup

	<b>Blockcycler</b>	<b>Rotor-Gene Q</b>
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 50-65°C*	10 sec, 95°C 15 sec, 50-65°C*
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase  Reporter Dye 1 (522 nm): FAM Reporter Dye 2 (553 nm): VIC (internal amplification control) Quencher Dye: TAMRA or BHQ  Passive Reference: none	Detection: End of extension phase Reporter Dye 1: Green Reporter Dye 2: Yellow (internal amplification control)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads">http://www.congen.de/en/company/downloads</a>		

\*Interactions of the oligonucleotides from the detection system and the internal amplification control are higher in the case of low annealing temperatures. We recommend using annealing temperatures of about 60°C.

3. Preparation of the PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the PCR instrument and start the run according to the setup.

**Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows amplifications in the detection system (FAM).

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system. In case of a **negative** result the inhibition control of the sample must be **positive**. Is this not the case the sample contains PCR-inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample needs to be improved.

**Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).