CONGEN

SureCC Color Compensation Kit IV

Art. No. F4012
3 calibration runs

User Manual

December 2025

Inhalt

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Dezember 2025

Content

1	General Information.....	11
1.1	Description.....	11
1.2	Kit components and storage	11
1.3	Additionally required equipment and materials	11
1.4	Precautions for users	11
2	Protocol	12
2.1	Preparation of the color compensation plate	12
2.2	Setup	13
2.2.1	program the detection format.....	13
2.2.2	Set a new experiment	14
2.2.3	Program the PCR layout – Subset Editor	15
2.2.4	Program the PCR layout – Sample Editor	16
2.3	Start of the experiment	16
2.4	Evaluation and creation of a CC-File.....	17
3	Further Information	17
3.1	Product Information.....	17
3.2	Technical Support.....	17

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1 General Information

1.1 Description

This kit is designed for generating color compensation objects specifically for 5plex real-time experiments on the Roche LightCycler® 480 II instrument. Using the five included solutions, users can perform a color compensation experiment. The resulting fluorescence data are processed to create a color compensation file, which corrects crosstalk between detection channels during 5plex real-time PCR runs. The generated color compensation object can thereafter be applied to analyze 5plex real-time PCR experiments on the Roche LightCycler® 480 II.

1.2 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Blank	1 x 400 µl	White
2	Dye 1	1 x 400 µl	Blue
3	Dye 2	1 x 400 µl	Green
4	Dye 3	1 x 400 µl	Yellow
5	Dye 4	1 x 400 µl	Orange
6	Dye 5	1 x 400 µl	Red

Store all reagents at –28 to -16°C and protected from light.

1.3 Additionally required equipment and materials

- real-time PCR instrument: LightCycler® 480 II
- real-time PCR consumable (plates, foils)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.4 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

2 Protocol

2.1 Preparation of the color compensation plate

For a color compensation experiment it is necessary to pipette 5 reactions with 20 µl for each dye and 5 reactions with 20 µl for the fluorescence background (blank). See Figure 1:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D	Blank		Dye 1		Dye 2		Dye 3		Dye 4		Dye 5	
E												
F												
G												
H												

Figure 1: pipetting scheme

Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Kit Code	Solution	Amount per reaction	Pipetting of 20 µl in following wells
1	Blank	20.0 µl	B1, C1, D1, E1, F1
2	Dye 1	20.0 µl	B3, C3, D3, E3, F3
3	Dye 2	20.0 µl	B5, C5, D5, E5, F5
4	Dye 3	20.0 µl	B7, C7, D7, E7, F7
5	Dye 4	20.0 µl	B9, C9, D9, E9, F9
6	Dye 5	20.0 µl	B11, C11, D11, E11, F11

Close the plate with an optical sealing foil after pipetting the solutions.

2.2 Setup

2.2.1 Program the detection format

Open the LightCycler® 480 II software. Please program under the icon for “tools” the necessary detection format (see Figure 2):


Figure 3

In the Tool window select “Detection Formats”. Click the “New” button to create a new detection format (see table below). Save it, for example as „Splex 440-488 | 465-510 | 533-580 | 533-610 | 618-660 “(see Figure 3).

Dye name	Filter Combination	Roche LightCycler® 480 II
ATTO		440-488
FAM		465 - 510
VIC		533 - 580
ROX		533 - 610
Cy5		618 - 660

Set the value for Quant- and Melt Factors as well as the Integration Time to 1 (see Figure 3).

Choose the “Close” button to exit the Tools window.

SureCC Color Compensation Kit IV (3 calibration runs)

Art. No. F4012

December 2025

2.2.2 Set a new experiment

After the programming of the detection format please choose the button „**New Experiment**“. The following window will open (see Figure 4):


Figure 4

As the first step select the correct detection format (please choose the correct name – see procedure in figures 2 and 3) and correct reaction volume (20 µl).

The following protocol steps are to program (see Figure 5):

(Please pay attention to the correct Cycle number und Analysis Mode.)


Figure 5

Please program following protocol steps, temperature targets and measurement points:

Program	Cycles / Analysis Mode	Temperature targets			
		Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp rate [°C/s]
Initial Denat.	1 / none	95	none	00:00:30	4.4
Cycling	5 / Quantification	95	none	00:00:15	4.4
		60	single	00:00:30	2.2
TM-Analyse	1 / Melting Curves	95	none	00:00:01	4.4
		50	none	00:00:30	2.2
		70	continuous		0.14 (Acquisitions per °C = 1)
Cooling	1 / none	50	none	00:00:01	2.2

The final screen should look like the screenshot below (see Figure 6):


Figure 6

2.2.3 Program the PCR layout – Subset Editor

In order to select the correct PCR layout please click on the **“Subset Editor”** button on the left side.

Please create a subset through a click on the **“Plus”**-button and give the subset a name for example **“CC Layout”**. Now please press the Strg-key and the left mouse button simultaneously. Mark all wells with the solution inside the plate. Then click on the **“Apply”** button and the subset is ready.

The final screen should look like the screenshot below (see Figure 7):


Figure 7

SureCC Color Compensation Kit IV

(3 calibration runs)

Art. No. F4012

December 2025

2.2.4 Program the PCR layout – Sample Editor

In the next step please select the “**Sample Editor**” button on the left side. Here please choose in the Step 1 **Select Workflow** the “**Color Comp**” tab. Now please select in Step 2 “**Select Samples**” dialog box the subset of the programming step before, e.g. “CC”. In the Dominant Channel dialog fields select the corresponding dyes (ATTO, FAM, VIC, ROX, CY5) in the different wells. For the reaction with the dye background (Blank), select “Water”. It is not necessary to give a “Sample Name” for the Dye’s (see Figure 8).


Figure 8

Reagent	Dominant Channel	Color
Blank	Water	Blank
Dye 1	440 - 488	ATTO
Dye 2	465 - 510	FAM
Dye 3	533 - 580	VIC
Dye 4	533 - 610	ROX
Dye 5	618 - 660	CY5

Please save your changes and the experiment by clicking on the disc icon (see Figure 8).

2.3 Start of the experiment

Place the PCR plate, which contains the reactions prepared before, into the LightCycler® 480 II.

Click the “**Start run**” button in the **Experiment** window (on the left side) to start the experiment.

2.4 Evaluation and creation of a CC-File

After the experiment ends, click the **"Analysis"** button in the left side of the window to open the **"Create New Analysis"** dialog box. Select **"Color Compensation"** and confirm your choice.

In the next dialog-box please select the subset und confirm your selection. Click now the button **"Calculate"** in the analysis window to perform color compensation analysis. Click **"Save CC Object"** button to save the Color Compensation file in the "CC" folder. The Color Compensation file is now ready to use.

To apply a Color Compensation file for data analysis, open under **"Analyse"** the respective multiplex assay and choose in **"Color Comb"** the needed filter combination. In the **"Color Comp"** drop-down menu (see figure 9) select "in Database" and choose the stored Color Compensation file you want to apply to the assay. The "Color Comp" button switches to **"Color Comp (On)"** to confirm that a Color Compensation is applied. The multiplex run can now evaluated.


Figure 9

Note: As the Color Compensation is instrument-specific, it is necessary to generate a Color Compensation object for every LightCycler® 480 II. A new object also has to be created after the optical system has been repaired.

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.