FAQ

SureFood® ALLERGEN

November 2024

CONGEN

 **Inhalt**

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November 2024

Content

1 Legal Regulation of Allergen Labeling..... 6

2 Which swabs are suitable for allergen analysis? 6

3 What must be considered with highly swelling samples? 6

4 Can all matrices be qualitatively analyzed for allergens using 4plex PCR? 6

5 Can all matrices be quantitatively analyzed for allergens using PCR? 6

6 Can the SureFood® Allergen QUANTARD 40 be used for spike experiments? 6

7 Why does the SureFood® Allergen QUANTARD 40 also contain milk and egg? 7

8 Does the DNA content of the sample or the SureFood® Allergen QUANTARD 40 need to be measured or adjusted? 7

9 How can the detection limit of the SureFood® ALLERGEN PCR Kits (Art. No. S3601 to S3620) be verified? 7

10 Can I dilute the QUANTARD 40 DNA to verify the detection limit of the SureFood® Allergen PCR Kits? 7

11 What is the measurement uncertainty of quantitative allergen PCR? 7

12 Is the Cp value the same for all allergens in the SureFood® Allergen QUANTARD 40? 7

13 What Cp fluctuations are possible between different SureFood® Allergen QUANTARD 40 lots? 7

14 What should be considered when evaluating the internal amplification control (IAC)? 7

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1 Legal Regulation of Allergen Labeling

According to current legislation, labeling is mandatory only for ingredients. Since 2014, the Food Information Regulation (EU) No. 1169/2011 governs the labeling of food allergens.

Manufacturers can label unintended contaminations, known as "cross-contacts." However, there is no legal regulation on the level at which this must occur. The labeling of traces entered as contamination is mainly for liability protection. For example, if an allergic person eats a muesli bar contaminated with peanuts and suffers an anaphylactic shock, the manufacturer could be held accountable if they knew about the peanut contamination and did not label it. This falls under the duty of care.

In terms of preventive consumer protection, official monitoring in Germany has established assessment values for allergens to uniformly evaluate analysis results. These values are based on the current scientific data of the VITAL 3.0 concept but are not legally binding. If the values exceed the assessment limits, it is checked whether it is an ingredient and, if necessary, the diligence of the company is reviewed.

2 Which swabs are suitable for allergen analysis?

For sampling in allergen analysis, we recommend sterile, dry viscose or cotton swabs (without a wooden stick). The swabs should be moistened with distilled water shortly before sampling to swab the relevant area. Then place the swab back in the tube. In PCR allergen analysis, no transport medium is needed. The DNA remains stable, and for DNA extraction, only the swab is cut off and lysed.

3 What must be considered with highly swelling samples?

For highly swelling samples, it may be necessary to add additional lysis buffer during DNA extraction. If SureFood® Allergen QUANTARD 40 is being extracted simultaneously for quantification, it is extracted as usual. The volume of the lysis buffer is not increased.

4 Can all matrices be qualitatively analyzed for allergens using 4plex PCR?

In mixed samples, there may be a loss of sensitivity in the detection channel with the lower concentration due to unequal mixing ratios (e.g., 99.9% almond and 0.1% cashew), especially if a Cp value before 20 is reached in one channel.

Therefore, the verification point "Limit of detection (LoD) in DNA mixtures" was included in the Validation Report. This is to ensure that a low-concentration analyte can be detected despite the background matrix. It should be noted that the background matrix must always be part of the overall system and form its own detection.

5 Can all matrices be quantitatively analyzed for allergens using PCR?

No, swabs, sponges, washing and cleaning solutions, and solid fats cannot be quantified for allergens using PCR.

Reason:

- Allergen quantification is performed with the external reference material SureFood® Allergen QUANTARD 40 and SureFood® PREP Advanced. This reference material is based on cornmeal and contains 40 mg of allergenic food component per kg of cornmeal. For analysis with external reference material, it is crucial that all sample preparation steps are carried out identically. Otherwise, the sample results are not comparable with those of the reference material. Therefore, only solid samples are quantified.
- Samples with high fat content require special treatment during DNA extraction. Solid samples with high fat content are not quantitatively analyzed because the SureFood® PREP Advanced is spin-filter-based, and fats can clog the spin filter. The fat content must first be removed, e.g., by hexane precipitation, which eliminates the requirement for identical sample preparation.
- It is also important that the samples are homogeneous. Otherwise, the quantitative results are not reproducible and meaningful.

6 Can the SureFood® Allergen QUANTARD 40 be used for spike experiments?

Yes, the SureFood® Allergen QUANTARD 40 can be used for conducting spiking experiments and internal verification.

7 Why does the SureFood® Allergen QUANTARD 40 also contain milk and egg?

The SureFood® Allergen QUANTARD 40 is a versatile external reference material. It was developed to accurately quantify allergenic components in foods from 1 mg to 400 mg per kg.

It contains all allergenic food components that must be labeled according to Regulation (EU) No. 1169/2011, except lactose and sulfite, at a concentration of 40 mg/kg.

8 Does the DNA content of the sample or the SureFood® Allergen QUANTARD 40 need to be measured or adjusted?

For qualitative and quantitative allergen analysis with the SureFood® ALLERGEN PCR Kits, it is not necessary to measure or adjust the DNA content of the sample or the reference material SureFood® ALLERGEN QUANTARD 40.

The total DNA content does not provide information about the proportion of food allergen DNA in a sample. The DNA content of SureFood® ALLERGEN QUANTARD 40 consists of a very small proportion of food allergen DNA and a significantly larger proportion of corn DNA. Therefore, measuring the DNA content is not helpful.

9 How can the detection limit of the SureFood® ALLERGEN PCR Kits (Art. No. S3601 to S3620) be verified?

For verifying the detection and determination limit, the SureFood® QUANTARD Allergen 1 ppm can be used as external reference material (Art. No. S3305).

10 Can I dilute the QUANTARD 40 DNA to verify the detection limit of the SureFood® Allergen PCR Kits?

The SureFood® Allergen QUANTARD 40 is an external reference material for the absolute quantification of food allergens. It consists of allergen-free cornmeal to which homogenized allergenic food components have been added. The goal was to create a reference material with a high background DNA content that closely resembles a real food matrix. Diluting the QUANTARD 40 DNA is not recommended, as it dilutes both the food allergen and the background matrix.

To verify the detection and determination limit of the SureFood® ALLERGEN Kits (Art. No. S3601 to S3620), the SureFood® QUANTARD Allergen 1 ppm (Art. No. S3305) is recommended.

11 What is the measurement uncertainty of quantitative allergen PCR?

The quantitative allergen PCR with SureFood® ALLERGEN Kits (Art. No. S3601 to S3620) has a measurement uncertainty of up to 30% (see validation reports). Near the determination limit, the measurement uncertainty can reach up to 100%.

12 Is the Cp value the same for all allergens in the SureFood® Allergen QUANTARD 40?

No, the Cp value of the SureFood® Allergen QUANTARD 40 varies from allergen to allergen. The reasons for this origin in the choice of target.

13 What Cp fluctuations are possible between different SureFood® Allergen QUANTARD 40 lots?

Each new batch of the SureFood® Allergen QUANTARD 40 is tested for homogeneity according to IUPAC.

During the production of the SureFood® ALLERGEN PCR Kits, the Cp values of the positive control, the standard DNA, and the SureFood® Allergen QUANTARD 40 are documented and monitored. For example, the Cp value for mustard in the SureFood® Allergen QUANTARD 40 remains constant. If the same DNA extraction is used for the sample and the reference material and the criteria for the standard curve from the manual are met, quantification is easily possible.

The determination in ppm or mg/kg is only possible in the range of 1-400 ppm mustard. Below the determination limit, a quantitative statement is not possible. In this case, the assessment would be positive below the determination limit.

14 What should be considered when evaluating the internal amplification control (IAC)?

The SureFood® PCR Kits have an internal amplification control (IAC). If inhibitory substances are present in the DNA, the signal of the amplification control is disturbed or the amplification is suppressed.

High concentrations of the amplicon, i.e., strongly positive samples, can also lead to a weak or absent signal of the IAC. A Cp value for the IAC is not necessary to obtain a positive result for the positive control.

The amplification of the IAC must only be evaluated for negative samples, i.e., when no amplification is visible in the detection channels.