
CONGEN

SureFast® PREP Bacteria

Art. No. F1021
100 extractions

User Manual

Efficient DNA preparation of bacteria

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Content

1	General Information	9
1.1	Description	9
1.2	Kit components and storage	9
1.3	Additionally required equipment and materials	10
1.4	Precautions for users	10
2	Protocol	11
2.1	Principle	11
2.2	Preparations	11
2.3	Extraction protocol	11
	1. Preparation of the basic material	11
	2. Lysis of the basic material	11
	3. Setting of optimal binding conditions	12
	4. Binding of the nucleic acids on a Spin Filter	12
	5. Purification of the bound nucleic acids	12
	6. Drying of the Spin Filter	12
	7. Elution of nucleic acids from the Spin Filter	12
3	Further Information	12
3.1	Product Information	12
3.2	Technical Support	12
4	Safety Information	13

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Kit Code	Reagenz/Material	Menge (je Box)
L	Lysis Buffer	1 x 25 ml
B	Binding Buffer*	1 x 4 ml
W	Wash Buffer**	1 x 18 ml
E	Elution Buffer	1 x 10 ml
S	Spin Filter	1 x 50 Filter
R	Receiver Tubes 2,0 ml	1 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes

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1 General Information

1.1 Description

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

The SureFast® PREP Bacteria kit has been validated and certified in combination with the SureFast® Salmonella PLUS real-time PCR kit (Art. No. F5111; AOAC PTM 041103), SureFast® Listeria 3plex ONE (Art. No. F5217; AOAC PTM 062501, MicroVal licence no. 2023LR114) and SureFast® STEC 4plex ONE, SureFast® Escherichia coli Serotype I 4plex, SureFast® Escherichia coli Serotype II 4plex (Art. No. F5265, F5167, F5168; AOAC PTM 052602).

In samples with a very low microbial activity (e.g. heat pasteurized egg white powder) the detection of viable or not culturable bacteria cells is possible. In case of other conditions the samples show a fast metabolization of these cells by accompanying flora.

To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing.

1.2 Kit components and storage

Kit Code	Reagent / Material	Amount (per box)
L	Lysis Buffer	1 x 25 ml
B	Binding Buffer*	1 x 4 ml
W	Wash Buffer**	1 x 18 ml
E	Elution Buffer	1 x 10 ml
S	Spin Filter	1 x 50 Filter
R	Receiver Tubes 2.0 ml	1 x 50 Tubes
T	Receiver Tubes 1.5 ml	1 x 50 Tubes

* Add isopropanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

** Add ethanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

The components are equally divided into two boxes.

All reagents of the kit should be stored dry and at room temperature (14-25°C).

1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase - 1.5 ml or 2.0 ml, with lock cap
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips, for volumes from 0.7 µL – 1000 µL
- Vortex mixer
- Thermomixer/ heating block (up to 99°C)
- micro centrifuge (up to 12,000 rpm) with a rotor for the reaction tubes
- ethanol for preparation of Wash Buffer (puriss, purity ≥96%)
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipuran® >99.7%; Applichem 2-Propanol for microbiology ; Sigma-Aldrich 2-Propanol)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

1.4 Precautions for users

- This test must only be performed by qualified and authorized laboratory personnel.
- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Strictly follow the working instructions.
- When handling samples, wear suitable protective equipment (disposable gloves, lab coat if necessary). After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

For more details see Material Safety Data Sheet, at www.congen.de/en/eifu/.

2 Protocol

2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

2.2 Preparations

General

Add 11 ml isopropanol to the Binding Buffer (**Code B**) and mix thoroughly. After adding isopropanol, the shelf life of the Binding Buffer is as indicated on the label, if stored at 14-25 °C.

Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly. After adding ethanol, the shelf life of the Wash Buffer is as indicated on the label, if stored at 14-25 °C.

Before each preparation

Preheating the Elution Buffer (**Code E**): Transfer the needed amount of Elution Buffer (up to 10% excess) into a reaction tube (not provided with the kit).

Pre-heat the Elution Buffer in a heating block to 60°C. The Elution Buffer is necessary for step 7.

2.3 Extraction protocol

1. Preparation of the basic material

Mix enrichment. Allow sample enrichment to settle for 5-10 minutes.

Transfer 1.0 ml of the enrichment or avulsions under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).

Centrifuge tubes for 5 min at 12,000 rpm.

Discard the fluid by aspirating with a micropipette and inactivate the supernatant per autoclaving.

For swabs, place them into a 1.5 ml reaction tube (not provided with the kit) and cut the shaft, so that you can close the cap of the extraction tube.

2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) to the reaction tube and cap tube. Mix it briefly by vortexing.

Place the tubes on a thermomixer with heating block and incubate under continuously shaking for 10 min at 99°C.

3. Setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Transfer 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

Note: Depending on the nature of the sample a smaller amount of liquid may be obtained after centrifugation. In this case take the actual volume for the following steps or increase the volume of the Lysis Buffer before lysis.

4. Binding of the nucleic acids on a Spin Filter

Add 200 µl Binding Buffer (**Code B**) to the supernatant and cap tube. Mix briefly by vortexing.

Place a Spin Filter (**Code S**) into a new 2.0 ml Receiver Tube (**Code R**). Transfer the solution onto the Spin Filter.

Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm.

After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

5. Purification of the bound nucleic acids

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

6. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature.

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

3 Further Information

3.1 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request
- Material Safety Data Sheet (Download: www.congen.de/en/eifu/)

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.

4 Safety Information

Lysis Buffer

Warning

H319-H412-P264-P273
P280-P305+P351+P338-
P337+P313-P501

- | | |
|------------------------|---|
| H319: | Causes serious eye irritation. |
| H412: | Harmful to aquatic life with long lasting effects. |
| P264: | Wash hands, forearms and face thoroughly after handling. |
| P273: | Avoid release to the environment. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P305+P351+P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313: | If eye irritation persists: Get medical advice/attention. |
| P501: | Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation. |

For further information we offer a Material Safety Data Sheet, see at www.congen.de/en/eifu/.
Alternatively please contact your distributor or send an e-mail to info@congen.de.